


FOOT & ANKLE Orthotics Catalog

Orthotic Solutions for
Management of Foot Drop &
Lower Extremity Instabilities



NaviGAIT™

For those requiring mild dorsiflexion assist

- Supports ankle dorsiflexion during swing phase
- Amount of dorsiflexion easily adjusted with the Velcro attachment on shin plate
- Low profile design makes it virtually invisible under slacks
- Plastic shin component can be trimmed with scissors

Item No.	Description	Height	Size
28410 0010	NaviGAIT	9"	X-Small
28410 0011	NaviGAIT	10.5"	Small/Medium
28410 0013	NaviGAIT	11.75"	Large/X-Large
28412 0010	Additional Strap w/Coil		X-Small
28412 0011	Additional Strap w/Coil		Small/Medium
28412 0013	Additional Strap w/Coil		Large/X-Large
28413 0010	Additional Strong Strap w/Coil		X-Small
28413 0011	Additional Strong Strap w/Coil		Small/Medium
28413 0013	Additional Strong Strap w/Coil		Large/X-Large
28414 0010	Extra Strap Set		X-Small
28414 0011	Extra Strap Set		Small/Medium
28414 0013	Extra Strap Set		Large/X-Large




4-Foot Wrap

Soft elastic forefoot wrap

- Makes it possible to wear NaviGAIT with:
 - Sandals
 - Shoes that do not have an anchor point
 - Barefoot

Item No.	Size	Circumference
28415 0011	Small	8.25" - 9"
28415 0012	Medium	9" - 10.25"
28415 0013	Large	10.25" - 11.25"





VISSION™ Stretch-Up

Nighttime treatment for plantar fasciitis

- Comfortable alternative to rigid night splints.
- Reduces the pain usually experienced during the first steps in the morning.
- Provides comfortable stretch to the calf and soft tissues of the foot
- Stretches the plantar fascia in an active and controlled way

Item No.	Size	Calf circumference
34108 0011	Small	< 7 ¾
34108 0012	Medium	7 ¾ - 15 ¾
34108 0013	Large	> 15 ¾



Night Splint Standard

Nighttime treatment to reduce the daytime pain associated with Plantar Fasciitis, Achilles Tendonitis, and other foot and ankle conditions.

- Gentle stretch for Achilles Tendon and Plantar Fascia
- Metatarsal wedge maximizes dorsiflexion stretch
- Adjustable, bilateral dorsiflexion straps
- Slip resistant sole
- Fits right or left foot

Item No.	Size	Shoe size Men	Shoe size Women
445730	Small	3 - 7	4 - 8
445740	Medium	8 - 10	9 - 11
445750	Large	11 +	12 +

M2 Ankle Laced Figure 8

Firm compressive support to help stabilize a weakened or injured ankle.

- Figure 8 strap replicates ankle taping
- Wide elastic cuff provides circumferential compression
- Oversized eyelets allow quick lace-up
- Design accommodates anatomical contours
- Cushioned achilles relief
- Convenient finger pulls on strap ends
- Fits right or left foot



Item No.	Size	Shoe size Men	Shoe size Women
445221	X-Small	3 - 5	4 - 6
445231	Small	6 - 7	7 - 9
445241	Medium	8 - 10	10 - 11
445251	Large	11 - 12	12 - 13
445261	X-Large	13 - 14	13 +

Ankle Foot Stabilizer

Stability and heat therapy for tendons

Ideal for patients with diabetes, arthritis, Raynaud's disease and other indications where warmth and stability is desired.

- Provides warmth so tendons are more responsive to therapy & exercise
- Includes medial and lateral stabilizers
- Adjustable closure strap for greater support of weak ankles
- Non-slip safety sole
- Fits right or left foot



Item No.	Size	Shoe Size Men	Shoe Size Women
47516 0011	Small	5 ½ - 7	6 ½ - 9
47516 0012	Medium	7 ½ - 10	9 ½ - 11
47516 0013	Large	10 ½ - 12	11 ½ - 13
47516 0014	X-Large	12 ½ - 15	13 ½ - 15




F.O.S. Foot Orthotic Shell

Optimize foot biomechanical function for mild to moderate pes planus/cavus

- High density EFM (Extra Firm Material)
- Microcell Puff(R) Interface is EVA foam, excellent shock absorber
- Grindable, Washable
- Trimmable and heat moldable
- Excellent shape retention and abrasion resistant

Item No.	Size	F.O.S Length	Adult Footplate Length
68000 0010	X-Small	8¾"	8¼"
68000 0011	Small	9½"	9"
68000 0012	Medium	10"	9⅝"
68000 0013	Large	10¾"	10⅝"
68000 0014	X-Large	11½"	11¼"
68000 0015	XX-Large	12"	N/A
			Pediatric Footplate Length
68372 0011	Small	6"	6 ¼"
68372 0012	Medium	6¾"	7 ⅛"
68372 0013	Large	7½"	7 ⅞"
68372 0014	X-Large	8¼"	N/A



WedgeKIT

Posting wedges for correcting or adjusting medial / lateral foot positioning

- Pressure sensitive Adhesive (PSA) Backed EVA Foam
- Includes one lateral wedge and 3 heel wedges
- Shore A60 Durometer

Color	Dimension	Slope
Grey	3.875"	10 mm - 1 mm
Chocolate	3.875"	7 mm - 1 mm
Black	3.875"	5 mm - 1 mm
Chocolate	7.75"	7 mm - 1 mm

Item No.	Description	Size
28766 0000	WedgeKIT	One Size

Vission Walker Clamshell Standard Tall

Circumferential compression without the weight

Our Vission Walker Clamshell Standard Tall is a unique clamshell design that offers circumferential compressive support and immobilization with a contoured-for-comfort liner.

- Rocker bottom with ample tread for more natural ambulation
- Extra heel width provides more surface contact and better stability
- Padded insole for better shock absorption
- Seamless foot liner for comfort
- Semi-rigid anterior and posterior panels accommodate anatomy
- Breathable, moisture wicking, contoured-for-comfort liner
- Steel reinforced rocker bottom



Item No.	Size	Shoe size Men	Shoe size Women
440020	X-Small	2 - 4	3 - 5½
440030	Small	4½ - 7	6 - 8
440040	Medium	7½ - 10½	8½ - 11½
440050	Large	11 - 12½	12 - 13½
440060	X-Large	13 +	14 +

Vission Walker Clamshell Air Tall

Pneumatic circumferential compression without the weight

Our Vission Clamshell Air Tall integrates a single pneumatic bladder with our unique anterior/posterior clamshell design to offer circumferential support and immobilization with a contoured-for-comfort liner.

- Integrated single pneumatic bladder system
- Rocker bottom with ample tread for more natural ambulation
- Extra heel width provides more surface contact and better stability
- Padded insole for better shock absorption
- Seamless foot liner for comfort
- Semi-rigid posterior and anterior panels accommodate anatomy
- Breathable, moisture wicking, contoured-for-comfort liner
- Steel reinforced rocker bottom



Item No.	Size	Shoe size Men	Shoe size Women
440120	X-Small	2 - 4	3 - 5½
440130	Small	4½ - 7	6 - 8
440140	Medium	7½ - 10½	8½ - 11½
440150	Large	11 - 12½	12 - 13½
440160	X-Large	13 +	14 +



Walker Standard Tall and Short

Compression and stability

Our Walker Standard Tall and Short offer circumferential compressive support and immobilization with a contoured-for-comfort liner.

- Rocker bottom with ample tread for more natural ambulation
- Extra heel width provides more surface contact and better stability
- Padded insole for better shock absorption
- Seamless foot liner for comfort
- Breathable, moisture wicking, contoured-for-comfort liner
- Steel reinforced rocker bottom

TALL

Item No.	Size	Shoe size Men	Shoe size Women
44032I	X-Small	2 - 4	3 - 5½
44033I	Small	4½ - 7	6 - 8
44034I	Medium	7½ - 10½	8½ - 11½
44035I	Large	11 - 12½	12 - 13½
44036I	X-Large	13 +	14 +



SHORT

Item No.	Size	Shoe size Men	Shoe size Women
44052I	X-Small	2 - 4	3 - 5½
44053I	Small	4½ - 7	6 - 8
44054I	Medium	7½ - 10½	8½ - 11½
44055I	Large	11 - 12½	12 - 13½
44056I	X-Large	13 +	14 +

Walker Air Tall and Short

Integrated pneumatic air system

Our Walker Air Tall and Short offer an integrated pneumatic single bladder air system combined with circumferential compressive support and immobilization with a contoured-for-comfort liner.

- Integrated single pneumatic bladder
- Rocker bottom with ample tread for more natural ambulation
- Extra heel width provides more surface contact and better stability
- Padded insole for better shock absorption
- Seamless foot liner for comfort
- Breathable, moisture wicking, contoured-for-comfort liner
- Steel reinforced rocker bottom



TALL

Item No.	Size	Shoe size Men	Shoe size Women
440421	X-Small	2 - 4	3 - 5½
440431	Small	4½ - 7	6 - 8
440441	Medium	7½ - 10½	8½ - 11½
440451	Large	11 - 12½	12 - 13½
440461	X-Large	13 +	14 +

SHORT

Item No.	Size	Shoe size Men	Shoe size Women
440621	X-Small	2 - 4	3 - 5½
440631	Small	4½ - 7	6 - 8
440641	Medium	7½ - 10½	8½ - 11½
440651	Large	11 - 12½	12 - 13½
440661	X-Large	13 +	14 +



Be aware that the brace may alter the patient's abilities for some activities, i.e. driving a car. Review the Instructions For Use supplied with the product with patient and provide them with a copy. Point out where discussed in IFU and remind the patient the importance of routine checking of orthosis for any signs of damage.

Product Warranty

Allard USA obligation under this warranty, is limited to the replacement or repair of any defective part or parts of the product within 30 days of fit. Normal wear and tear during the usage of a product is not considered a defect. It does not cover damage by tools, application of heat, or excessive bending/flexing of brace.

Patient Satisfaction Guarantee

If the wearer is not totally satisfied with the product initially fitted within 10 days, Allard USA will at the provider's discretion, provide a FREE alternative product for the wearer to try, or issue full invoice cost credit.

Return Policy

If you have purchased the product from Allard USA:

To help us serve you better and assure you receive appropriate credit, please call toll free 888-678-6548 to receive a Return Authorization (RA) number prior to shipping the return.

If you have purchased the product from a distributor:

Please contact that distributor to obtain a Return Authorization.

If you have purchased the product from any company other than Allard USA:

Please contact that company to obtain a Return Authorization

Delivery Times

It is our goal to ship all orders received by 4:00pm, same day. If for any reason we are unable to meet that goal, you will be contacted.

Terms

There is no minimum order. There is a \$6.00 handling charge per order.

Payment terms are net 30 days. Overdue invoices are subject to a charge of 1.5% of the unpaid balance or .50% per month, whichever is greater.

Ordering Information

To place an order, please contact Customer Service:

Call toll free
888-678-6548

Fax
800-289-0809

Email
orders@allardusa.com

Address

Allard USA, Inc.
300 Forge Way, Suite 3
Rockaway, NJ 07866-2056

Hours of Operation

Monday through Friday
8:30AM - 6:00 PM EST

Support for Better Life!



allardusa.com

ALLARD USA, INC.
300 Forge Way, Suite 3
Rockaway, NJ 07866-2056

info@allardusa.com
Toll Free 888-678-6548
Fax 800-289-0809